


Reference: BS131.800

## SCANNER CD60

HPTLC scanning densitometer


- Absorbance or fluorescence measurement**
- Remission or transmission mode**
- For objects up to 265 x 200 x 4 mm**
- Automatic starting and switching of filters and lamps**
- Rapid data collection and evaluation**
- Recording spectra from 190 - 900 nm**
- Automatic background noise correction**
- Software-controlled by ProQuant**
- Ease of operation**
- Reproducible and reliable results**



For quantitative determination of samples, the HPTLC-Densitometer CD60 converts the spots/bands of the single substances into a chromatogram curve. It measures the absorbance or fluorescence of separated compounds in transmission or reflection mode. The HPTLC-Densitometer CD60 is controlled by ProQuant software which also enables quantitative evaluation of the generated data.

The HPTLC-Densitometer CD60 works within a spectral range of 190 - 900 nm. This is provided by three light sources: a deuterium lamp (190 to 340 nm), a halogen lamp (340 to 900 nm) as well as a mercury lamp. Once the wavelength is selected, the densitometer will automatically start to scan the entire plate. It measures the absorbance or fluorescence reflected or transmitted by each sample.

This will be stored in the software in the form of peak tables. These tables consist of R<sub>f</sub> values and area of each spot. Therefore, you can carry forward the quantitative evaluation of the generated densitometric data by ProQuant software



## MEASUREMENT AND EVALUATION

### Method Types

Method for chromatogram  
Method for multi-wavelength scan  
Method for spectrum

### Recording Modes


Remission and transmission  
Absorbance or fluorescence  
Linear and Meander scan  
Two-wavelength measurement  
Multi-wavelength measurement

### Results


Peak lists  
Results for sample and standard  
Automatic integration with manual correction facility  
Linear, polynomial or Michaelis-Menten function

## TECHNICAL SPECIFICATIONS


Object size: Up to 265 x 200 x 4 mm  
Spectral range: 190 - 900 nm  
Filters: 370, 420, 450, 550, orange, UV  
Max. scan length: 5 to 195 mm  
Max. scan width: 5 to 260 mm  
Slit width: 0.4 to 10 mm  
Slit height: 0.02 to 2 mm  
Dimensions: (W x H x D) 730 x 550 x 300 mm  
Weight: 30 kg



« Multi-wavelength » mode



« Chromatogram » mode



« Overlay » mode.

## REFERENCES FOR ORDER

Reference	Description
BS131.800	HPTLC-Densitometer CD60, 230 V, incl. interface box, software ProQuant
BS131.801	HPTLC-Densitometer CD60, 110 V, incl. interface box, software ProQuant
BS131.816	Software Provalid, program for automatic validation
BS131.830	Software Spectra Calc, program for compilation of spectra libraries
BS131.825	IQ/OQ documents for HPTLC-Densitometer CD60

Technical changes reserved. Bionis is not responsible for any misprints, errors which may result in any losses, claims or costs